

PRODUCT INFORMATION

Target	SNAPC3
Synonyms	PTFbeta; SNAP50
Description	Purified recombinant protein of Human small nuclear RNA activating complex, polypeptide 3, 50kDa (SNAPC3), full length, with N-terminal HIS tag, expressed in E. coli, 50ug
Delivery	1 week
Uniprot ID	Q92966
Expression Host	E. coli
Tag	N-His
Molecular Characterization	N/A
Molecular Weight	46.6 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25mM Tris, pH8.0, 150mM NaCl, 10% glycerol, 1% Sarkosyl.
Storage & Shipping	Store at -80°C.
Background	Part of the SNAPc complex required for the transcription of both RNA polymerase II and III small-nuclear RNA genes. Binds to the proximal sequence element (PSE), a non-TATA-box basal promoter element common to these 2 types of genes. Recruits TBP and BRF2 to the U6 snRNA TATA box.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

