

PRODUCT INFORMATION

Tag	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 C-Flag&Strep Tag
Target	TAS1R2
Synonyms	GPR71; T1R2; TR2
Description	Human TAS1R2-Strep full length protein-synthetic nanodisc
Delivery	6~8weeks
Uniprot ID	Q8TE23
Expression Host	HEK293
Protein Families	Druggable Genome, Transmembrane
Protein Pathways	Taste transduction
Molecular Weight	The human full length TAS1R2-Strep protein has a MW of 95.2 kDa
Formulation & Reconstitution	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Putative taste receptor. TAS1R2/TAS1R3 recognizes diverse natural and synthetic sweeteners.
Usage	Research use only
Conjugate	Unconjugated

