

PRODUCT INFORMATION

Target	TCL1B
Synonyms	SYN-1; TML1
Description	Recombinant protein of human T-cell leukemia/lymphoma 1B (TCL1B), transcript variant 1
Delivery	2-3 weeks
Uniprot ID	O95988
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	14.7 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Enhances the phosphorylation and activation of AKT1 and AKT2.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

DIMA BIO CONFIDENTIAL

