

PRODUCT INFORMATION

Target	TGFA
Synonyms	TGF-alpha;ETGF;TGF type 1
Description	Recombinant Human TGFA Protein with C-terminal human Fc tag
Delivery	In Stock
Uniprot ID	P01135
Expression Host	HEK293
Tag	C-Human Fc Tag
Molecular Characterization	TGFA(Glu24-Gln98) hFc(Glu99-Ala330)
Molecular Weight	The protein has a predicted molecular mass of 34.1 kDa after removal of the signal peptide. The apparent molecular mass of TGFA-hFc is approximately 35-55 kDa due to glycosylation.
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	This gene encodes a growth factor that is a ligand for the epidermal growth factor receptor, which activates a signaling pathway for cell proliferation, differentiation and development. This protein may act as either a transmembrane-bound ligand or a soluble ligand. This gene has been associated with many types of cancers, and it may also be involved in some cases of cleft lip/palate. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Human TGFA Protein, hFc Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

