

PRODUCT INFORMATION

Target	TLX3
Synonyms	HOX11L2; RNX
Description	Recombinant protein of human T-cell leukemia homeobox 3 (TLX3)
Delivery	2-3 weeks
Uniprot ID	O43711
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	31.7 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	The protein encoded by this gene is an orphan homeobox protein that encodes a DNA-binding nuclear transcription factor. A translocation [t(5;14)(q35;q32)] involving this gene is associated with T-cell acute lymphoblastic leukemia (T-ALL) in children and young adults. [provided by RefSeq, Nov 2015]
Usage	Research use only
Conjugate	Unconjugated

