

PRODUCT INFORMATION

|                              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Tag                          | <b>Warning:</b> Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806<br>C-Flag Tag                                                                                                                                                                                                                                                                                                                    |
| Target                       | TRPM6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Synonyms                     | CHAK2, HMGX, HOMG, HOMG1, HSH                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Description                  | Human TRPM6 full length protein-synthetic nanodisc                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Delivery                     | 6~8weeks                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Uniprot ID                   | Q9BX84                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Expression Host              | HEK293                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Protein Families             | Ion Channels: Transient receptor potential                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Protein Pathways             | N/A                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Molecular Weight             | The human full length TRPM6 protein has a MW of 231.7kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Formulation & Reconstitution | Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.                                                                                                                                                                                                                                                                                                  |
| Storage & Shipping           | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.                                                                                                                                                                                                                                                                                               |
| Background                   | This gene is predominantly expressed in the kidney and colon, and encodes a protein containing an ion channel domain and a protein kinase domain. It is crucial for magnesium homeostasis, and plays an essential role in epithelial magnesium transport and in the active magnesium absorption in the gut and kidney. Mutations in this gene are associated with hypomagnesemia with secondary hypocalcemia. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Apr 2016] |
| Usage                        | Research use only                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Conjugate                    | Unconjugated                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |

