

PRODUCT INFORMATION

Target	TSNAXIP1
Synonyms	TXI1
Description	Recombinant protein of human translin-associated factor X interacting protein 1 (TSNAXIP1)
Delivery	2-3 weeks
Uniprot ID	Q2TAA8
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	76.6 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	Possible role in spermatogenesis.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

