

PRODUCT INFORMATION

Target	TSPY3
Synonyms	CT78; TSPY; TSPY1; TSPY10
Description	Recombinant protein of human testis specific protein, Y-linked 3 (TSPY3)
Delivery	2-3 weeks
Uniprot ID	Q01534
Expression Host	HEK293T
Tag	C-Myc/DDK
Molecular Characterization	N/A
Molecular Weight	34.9 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Storage & Shipping	Store at -80°C.
Background	May be involved in sperm differentiation and proliferation.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

