

PRODUCT INFORMATION

Target	Z DNA binding protein
Synonyms	C20orf183; DAI; DLM-1; DLM1
Description	Purified recombinant protein of Human Z-DNA binding protein 1 (ZBP1), transcript variant 1, full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Delivery	1 week
Uniprot ID	Q9H171
Expression Host	E. coli
Tag	N-His
Molecular Characterization	N/A
Molecular Weight	46.2 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25mM Tris, pH8.0, 150 mM NaCl, 10% glycerol, 1 % Sarkosyl
Storage & Shipping	Store at -80°C.
Background	This gene encodes a Z-DNA binding protein. The encoded protein plays a role in the innate immune response by binding to foreign DNA and inducing type-I interferon production. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2011]
Usage	Research use only
Conjugate	Unconjugated

