

PRODUCT INFORMATION

Target	ZNF500
Synonyms	ZKSCAN18; ZSCAN50
Description	Purified recombinant protein of Human zinc finger protein 500 (ZNF500), full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
Delivery	2-3 weeks
Uniprot ID	O60304
Expression Host	E. coli
Tag	N-His
Molecular Characterization	N/A
Molecular Weight	53.5 kDa
Purity	> 80% as determined by SDS-PAGE and Coomassie blue staining
Formulation & Reconstitution	25mM Tris, pH8.0, 150 mM NaCl, 10% glycerol, 1 % Sarkosyl
Storage & Shipping	Store at -80°C.
Background	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated

