

**PRODUCT INFORMATION**

<b>Target</b>	ZNF669
<b>Synonyms</b>	N/A
<b>Description</b>	Purified recombinant protein of Human zinc finger protein 669 (ZNF669), transcript variant 1, full length, with N-terminal HIS tag, expressed in E.Coli, 50ug
<b>Delivery</b>	2-3 weeks
<b>Uniprot ID</b>	Q96BR6
<b>Expression Host</b>	E. coli
<b>Tag</b>	N-His
<b>Molecular Characterization</b>	N/A
<b>Molecular Weight</b>	52.4 kDa
<b>Purity</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Formulation &amp; Reconstitution</b>	50mM Tris, 8M Urea, pH8.0.
<b>Storage &amp; Shipping</b>	Store at -80°C.
<b>Background</b>	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated

