

**PRODUCT INFORMATION**

<b>Target</b>	ADAM9
<b>Synonyms</b>	AU020942;MDC9;mKIAA0021;Mltng
<b>Description</b>	Recombinant mouse ADAM9 protein with C-terminal 6×His tag
<b>Delivery</b>	In Stock
<b>Uniprot ID</b>	Q61072
<b>Expression Host</b>	HEK293
<b>Tag</b>	C-6×His Tag
<b>Molecular Characterization</b>	Mouse ADAM9(Gly30-Asp697) 6×His tag
<b>Molecular Weight</b>	The protein has a predicted molecular mass of 74.3 kDa after removal of the signal peptide. The apparent molecular mass of mADAM9-His is approximately 55-70 kDa due to glycosylation.
<b>Purity</b>	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Yefei_Storage</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	Cleaves and releases a number of molecules with important roles in tumorigenesis and angiogenesis, such as TEK, KDR, EPHB4, CD40, VCAM1 and CDH5 (PubMed:19273593). May mediate cell-cell, cell-matrix interactions and regulate the motility of cells via interactions with integrins (PubMed:10825303).[UniProtKB/Swiss-Prot Function]
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated






Figure 1. Mouse ADAM9 Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

