

PRODUCT INFORMATION

Target	CD48
Synonyms	BCM1;SLAMF2;BLAST;BLAST1;MEM-102;TCT.1;BCM-1;SLAMF-2;BLAST-1
Description	Recombinant mouse CD48 protein with C-terminal human Fc tag
Delivery	In Stock
Uniprot ID	P18181
Expression Host	HEK293
Tag	C-Human Fc Tag
Molecular Characterization	Mouse CD48(Phe23-Ser217) hFc(Glu99-Ala330)
Molecular Weight	The protein has a predicted molecular mass of 48.5 kDa after removal of the signal peptide. The apparent molecular mass of mCD48-hFc is approximately 55-70 kDa due to glycosylation.
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Ligand for CD2. Might facilitate interaction between activated lymphocytes. Probably involved in regulating T-cell activation.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated




Figure 1. Mouse CD48 Protein, hFc Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

