

**PRODUCT INFORMATION**

<b>Target</b>	CDH17
<b>Synonyms</b>	HPT-1;HPT-1/LI
<b>Description</b>	Recombinant mouse CDH17 protein with C-terminal 6×His tag
<b>Delivery</b>	In Stock
<b>Uniprot ID</b>	Q9R100
<b>Expression Host</b>	HEK293
<b>Tag</b>	C-6×His Tag
<b>Molecular Characterization</b>	Mouse CDH17(Phe26-Met786) 6×His tag
<b>Molecular Weight</b>	The protein has a predicted molecular mass of 85.1 kDa after removal of the signal peptide. The apparent molecular mass of mCDH17-His is approximately 70-130 kDa due to glycosylation.
<b>Purity</b>	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Yefei_Storage</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	Enables proton-dependent oligopeptide secondary active transmembrane transporter activity. Involved in germinal center B cell differentiation; marginal zone B cell differentiation; and spleen development. Acts upstream of or within calcium-dependent cell-cell adhesion via plasma membrane cell adhesion molecules; homophilic cell adhesion via plasma membrane adhesion molecules; and oligopeptide transport. Located in cell surface and plasma membrane. Is expressed in several structures, including common bile duct; gut; spleen; trunk; and urogenital sinus. Orthologous to human CDH17 (cadherin 17). [provided by Alliance of Genome Resources, Apr 2022]
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated






Figure 1. Mouse CDH17 Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

