

PRODUCT INFORMATION

Target	TSLP
Synonyms	TSLP
Description	Recombinant mouse TSLP protein with C-terminal human Fc tag
Delivery	In Stock
Uniprot ID	Q9JIE6
Expression Host	HEK293
Tag	C-Human Fc Tag
Molecular Characterization	Mouse TSLP(Tyr20-Glu140) hFc(Glu99-Ala330)
Molecular Weight	The protein has a predicted molecular mass of 40.1 kDa after removal of the signal peptide. The apparent molecular mass of mTSLP-hFc is approximately 35-55 kDa due to glycosylation.
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Yefei_Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Enables cytokine activity. Predicted to be involved in several processes, including defense response to other organism; positive regulation of cytokine production; and positive regulation of signal transduction. Predicted to be located in extracellular region. Is expressed in ductus deferens; epididymis; ileum; and testis. Used to study atopic dermatitis; autoimmune hemolytic anemia; cryoglobulinemia; and membranoproliferative glomerulonephritis. Orthologous to human TSLP (thymic stromal lymphopoietin). [provided by Alliance of Genome Resources, Apr 2022]
Usage	Research use only
Conjugate	Unconjugated





Figure 1. Mouse TSLP Protein, hFc Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

