

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 H8, mAb5T4
Conjugate	PE-conjugated
Synonyms	M6P1;WAIF1;TPBG
Applications	Flow Cyt
Recommended Dilutions	1µl/test
Formulation & Reconstitution	Liquid;PBS with,0.18%BSA;0.1%Proclin300
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	5T4
Uniprot ID	Q13641
Description	PE-conjugated Anti-5T4 (H8 biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only




Figure 1. Flow cytometry analysis with 1 µl/test PE-conjugated Anti-5T4 (H8 biosimilar) mAb on HEK293 cells (Blue histogram) or Jurkat cells(Red histogram).

