Cat. No. DMC101149P



PRODUCT INFORMATION

Warning . Undefined variable \$hasAttributeValueDescription in Calwawroottmirror.dimablo.com/wp-content\plugins/woocommerce-print-products-public.php on line 2806 2612 Clone ID

APN;CD13;GP150;LAP1;P150;PEPN Synonyme

Host Species Rabbit

PE-conjugated Anti-ANPEP antibody(26H2), IgG1 Chimeric mAb Description

Delivery Under Development

Uniprot ID

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Flow Cyt Applications Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Liquid PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping

Store at 2°C-8°C for 6 months

Email: info@dimabio.com Website: www.dimabio.com

Background

Conjugate

PE-conjugated
All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

