

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806
Conjugate	PE-conjugated
Synonyms	BAFFR;TNFRSF13C;BAFF-R;BROMIX;CD268;CVID4;prolixin
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1µl/test
Formulation & Reconstitution	Liquid,PBS with,0.18%BSA,0.1%Proclin300
Host Species	Homo sapiens
IgG type	IgG1
Reactivity	Human
Target	BAFF-R
Uniprot ID	Q96RJ3
Description	PE-conjugated Anti-BAFFR (ianalumab biosimilar) mAb
Delivery	In Stock
Yefei_Storage	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only




Figure 1. Flow cytometry analysis with 1 µl/test PE-conjugated Anti-BAFFR(ianalumab biosimilar) on HEK293 cells transfected with Human BAFFR protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

