

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 BlueBird HuC11D5.3 (bb2121)
Conjugate	PE-conjugated
Synonyms	CD269;TNFRSF17;BCM;BCMA
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:250
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	BCMA
Uniprot ID	Q02223
Description	PE-conjugated Anti-BCMA (biosimilar) mAb
Delivery	In Stock
Yefei_Storage	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only




Figure 1. Flow cytometry analysis with 100 µl/test (1:250) PE-conjugated Anti-BCMA (biosimilar) mAb (BME100016P) (Blue histogram) or isotype control mAb (Red histogram) on K562-BCMA cells.

