Cat. No. DMC101034P



## **PRODUCT INFORMATION**

Warning: Undefined variable ShacktributeValueDescription in C:\text{Newworot\text{mirror.dimabio.com\wp-content\plugins\woocommerce-print-products-public.php on line 2806 302. Clone ID

BTN342

BT3.2;BTF4;BTN3.2;CD277 Synonyme

Host Species Rabbit

PE-conjugated Anti-BTN3A2 antibody(3G2); IgG1 Chimeric mAb Description

Delivery Under Development

Uniprot ID

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Flow Cyt Applications Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Liquid PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping

ay class 1

A sequencing or reverse This gene encodes a member of the immunoglobulin superfamily, which esides in the juxta-telomeric region of the major histocompatability class 1 Alternatively soliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq. Jun 2013] Research use only Background

Conjugate PE-conjugated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

Email: info@dimabio.com Website: www.dimabio.com



