

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 ICP-B04;ICP B04;MOR 208;MOR208;MOR-208;XENP-5574;XENP5574;XENP 5574;XmAb-5574;XmAb5574;XmAb 5574
Conjugate	PE-conjugated
Synonyms	CD19;B4;CVID3;MGCI2802
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1µl/test
Formulation & Reconstitution	Liquid,PBS with,0.18%BSA,0.1%Proclin300
Host Species	N/A
IgG type	IgG1
Reactivity	Human
Target	CD19
Uniprot ID	P15391
Description	PE-conjugated Anti-CD19(tafasitamab biosimilar) mAb
Delivery	In Stock
Yefei_Storage	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only




Figure 1. Flow cytometry analysis with 1µl/test PE-conjugated Anti-CD19(tafasitamab biosimilar) mAb on Raji cells (Blue histogram) or 293T cells (Red histogram).

