

## **PRODUCT INFORMATION**

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 bM211 Clone ID

Synonyme

Host Species Rabbit

PE-conjugated Anti-CD43 antibody(DM211); Rabbit mAb Description

Delivery Under Development Uniprot ID P16150 IgG type Rabbit IgG Clonality Monoclonal Reactivity Human Flow Cyt Applications

Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Liquid□PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping

Email: info@dimabio.com Website: www.dimabio.com

Background

Usage Research use only Conjugate

PEC-onjugated
All DIMA recombinant antibodies are genuinely generated by DIMA Blotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

