Cat. No. DMC100429P



PRODUCT INFORMATION

Warning: Undefined variable shasAttributeValueDescription in Cs\www.root\wirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public.chass-woocommerce-print-products-public.php on line 2806
bMC429 Clone ID

CD81

Synonyme

Host Species Rabbit

PE-conjugated Anti-CD81 antibody(DMC429); IgG1 Chimeric mAb Description

Delivery Under Development

Uniprot ID

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Flow Cyt Applications Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Liquid□PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping Store at 2°C-8°C for 6 months

Background

Conjugate PE-conjugated

protein seq. All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

> Email: info@dimabio.com Website: www.dimabio.com

