Cat. No. DMC100285P



PRODUCT INFORMATION

Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 bMC285 Clone ID

CD96 TACTILE Synonyme Host Species Rabbit

PE-conjugated Anti-CD96 antibody(DMC285); IgG1 Chimeric mAb Description

Delivery Under Development

Uniprot ID IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Flow Cyt Applications Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Liquid PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping Store at 2°C-8°C for 6 months

The protein encoded by this gene belongs to the immunoglobulin superfamily. It is a type I membrane protein. The protein may play a role in the adhesive interactions of activated T and NK cells during the late phase of the immune response. It may also function in antigen presentation. Alternative splicing generates multiple transcript variants encoding distinct isoforms.

Research use only Background

> Email: info@dimabio.com Website: www.dimabio.com

Conjugate PE-conjugated

...de in the ...door. Alternative. All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.



