

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 CDH17-ADC-653-H43
Conjugate	PE-conjugated
Synonyms	Cadherin-17
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1µl/test
Formulation & Reconstitution	Liquid,PBS with,0.18%BSA,0.1%Proclin300
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	CDH17
Uniprot ID	Q12864
Description	PE-conjugated Anti-CDH17(TORL-3-600 biosimilar) mAb
Delivery	Under Development
Yefei_Storage	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only




Figure 1. Flow cytometry analysis with 1 µl/test PE-conjugated Anti-CDH17(TORL-3-600 biosimilar) on SNU-5 cell line (Blue histogram) or 293T cell line (Red histogram).

