Cat. No. DMC100675P



## **PRODUCT INFORMATION**

Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 DMG/0765 Clone ID

Synonyme

Host Species Rabbit

PE-conjugated Anti-CDH6 antibody(DMC675); IgG1 Chimeric mAb Description

Delivery Under Development

Uniprot ID AAH00019

IgG type Rabbit/Human Fc chimeric IgG1 Clonality Monoclonal

Reactivity Human Flow Cyt Applications Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Liquid PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping

This gene encodes a member of the cadherin superfamily. Cadherins are membrane glycoproteins that mediate homophilic cell-cell adhesion and play critical roles in cell differentiation and morphogenesis. The encoded protein is a type II cadherin and may play a role in kidney development as well as endometrium and placenta formation. Decreased expression of this gene may be associated with tumor growth and metastasis. [provided by RefSeq: May 2011]

Research use only

Email: info@dimabio.com Website: www.dimabio.com

Usage Coniugate

PE-conjugated
All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

