

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 Chongqing Precision
Conjugate	PE-conjugated
Synonyms	CEACAM-5, CD66e, CEA, Meconium antigen 100
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	CEACAM5
Uniprot ID	P06731
Description	PE-conjugated Anti-CEACAM5 (Chongqing Precision biosimilar) mAb
Delivery	Under Development
Yefei_Storage	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMABIO CONFIDENTIAL

