

## PRODUCT INFORMATION

|   |  |
|---|--|
| <b>Common Name</b>                      | <b>Warning:</b> Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806<br>BMS-936564, MDX-1338 |
| <b>Conjugate</b>                        | PE-conjugated  |
| <b>Synonyms</b>                         | CD184;D2S201E;FB22;HM89;HSY3RR;LAP-3;LAP3;LCR1;LESTR;NPY3R;NPYR;NPYRL;NPYY3R;WHIM;WHIMS  |
| <b>Applications</b>                     | Flow Cyt   |
| <b>Recommended Dilutions</b>            | Flow Cyt 1:100   |
| <b>Formulation &amp; Reconstitution</b> | Liquid PBS with 0.05% Proclin300, 1% BSA   |
| <b>Host Species</b>                     | Homo sapiens   |
| <b>IgG type</b>                         | IgG4   |
| <b>Reactivity</b>                       | Human  |
| <b>Target</b>                           | CXCR4  |
| <b>Uniprot ID</b>                       | P61073   |
| <b>Description</b>                      | PE-conjugated Anti-CXCR4(ulocuplumab biosimilar) mAb   |
| <b>Delivery</b>                         | Under Development  |
| <b>Yefei_Storage</b>                    | Store at 2°C-8°C for 6 months  |
| <b>Background</b>                       | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.  |
| <b>Usage</b>                            | Research use only  |

DIMABIO CONFIDENTIAL

