

## PRODUCT INFORMATION

<b>Common Name</b>	<b>Warning:</b> Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 Fab C225,IMC-225
<b>Conjugate</b>	PE-conjugated
<b>Synonyms</b>	EGFR;ERBB;ERBB1;HER1;PIG61;mENA
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 0.5µl/test
<b>Formulation &amp; Reconstitution</b>	Liquid,PBS with,0.18%BSA,0.1%Proclin300
<b>Host Species</b>	Chimeric
<b>IgG type</b>	IgG1
<b>Reactivity</b>	Human
<b>Target</b>	EGFR
<b>Uniprot ID</b>	P00533
<b>Description</b>	PE-conjugated Anti-EGFR (Cetuximab biosimilar) mAb
<b>Delivery</b>	In Stock
<b>Yefei_Storage</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only




Figure 1. Flow cytometry analysis with 0.25 µl/test PE-conjugated Anti-EGFR (Cetuximab biosimilar) mAb on SNU-5 cells (Blue histogram) or Jurkat cells (Red histogram).

