

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\publicclass-woocommerce-print-products-public.php on line 2806 XENP28278
Conjugate	PE-conjugated
Synonyms	CD203c;NPP3;PD-1beta
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:400
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	ENPP3
Uniprot ID	O14638
Description	PE-conjugated Anti-ENPP3(Xencor) mAb
Delivery	In Stock
Yefei_Storage	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only




Figure 1. Flow cytometry analysis with 100 µl/test (1:400) PE-conjugated Anti-ENPP3(XENP28278) mAb (BME100151P) on HEK293 cells transfected with human ENPP3 (Blue histogram) or HEK293 cells(Red histogram).

