

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 BFCR-4350-A, BFCR-4350A, RG-6160, RO-7187797
Synonyms	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 BXMAS1;CD307;CD307e;FCRH5;IRTA2;PRO820
Conjugate	PE-conjugated
Applications	Flow Cyt
Delivery	In Stock
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	FCRL5
Uniprot ID	Q96RD9
Description	PE-conjugated Anti-FCRL5(cevostamab biosimilar without anti-CD3 portion) mAb
Recommended Dilutions	Flow Cyt 1:200
Yefel, Storage	Store at 2°C-8°C for 6 months
Formulation & Reconstitution	Liquid,PBS with 0.1% Proclin300, 0.18% BSA
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only




Figure 1. Flow cytometry analysis with 100 μ l/test (1:200) PE-conjugated Anti-FCRL5(cevostamab biosimilar without anti-CD3 portion) mAb (BME100089P) on HEK293 cells transfected with human FCRL5 (Blue histogram) or HEK293 cells (Red histogram).

