

PRODUCT INFORMATION

Clone ID	8D1
Target	GAS6
Synonyms	AXLLG;AXSF
Host Species	Rabbit
Description	PE-conjugated Anti-GAS6 antibody(8D1), IgG1 Chimeric mAb
Delivery	3-4 weeks
Uniprot ID	Q14393
IgG type	Rabbit/Human Fc chimeric IgG1
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	This gene encodes a gamma-carboxyglutamic acid (Gla)-containing protein thought to be involved in the stimulation of cell proliferation. This gene is frequently overexpressed in many cancers and has been implicated as an adverse prognostic marker. Elevated protein levels are additionally associated with a variety of disease states, including venous thromboembolic disease, systemic lupus erythematosus, chronic renal failure, and preeclampsia. [provided by RefSeq, Aug 2014]
Usage	Research use only
Conjugate	PE-conjugated

