PE-conjugated Anti-GPC1 antibody(1B8); IgG1 Chimeric mAb

Cat. No. DMC101004P



PRODUCT INFORMATION

Warning . Undefined variable \$hasAttributeValueDescription in Calwawroottmirror.dimablo.com/wp-content\plugins/woocommerce-print-products-public.php on line 2806
188 Clone ID

GPC1 Glypican 1;PC1 Synonyme Host Species Rabbit

PE-conjugated Anti-GPC1 antibody(1B8); IgG1 Chimeric mAb Description

Delivery Under Development

Uniprot ID

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Flow Cyt Applications Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Liquid PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping Store at 2°C-8°C for 6 months

Cell surface heparan sulfate proteoplycans are composed of a membrane-associated protein core substituted with a variable number of heparan sulfate chains. Members of the glypical-related integral membrane proteoglycan family (GRIPS) contain a core protein anchored to the cytoplasmic membrane via a glycosyl phosphatidylinositol linkage. These proteins may play a role in the control of cell division and growth regulation. [provided by RefSeq, Jul 2008]

Email: info@dimabio.com Website: www.dimabio.com

Research use only Usage Coniugate

PE-conjugated
All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

