

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 AMG888_U3-1287
Conjugate	PE-conjugated
Synonyms	HER3:ERBB3
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:400
Formulation & Reconstitution	Liquid PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG1
Reactivity	Human
Target	HER3
Uniprot ID	P21860
Description	PE-conjugated Anti-HER3 (patritumab biosimilar)mAb
Delivery	In Stock
Yefei_Storage	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only




Figure 1. Flow cytometry analysis with 100 µl/test (1:400) PE-conjugated Anti-HER3 (patritumab biosimilar)mAb (BME100057P) on HEK293 cells transfected with human HER3 (Blue histogram) or HEK293 cells (Red histogram).

