

PRODUCT INFORMATION

Common Name	RO4858696-000,teprotumumab-trbw
Synonyms	CD221, IGFIR, IGFR, JTK13
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Homo sapiens
IgG type	IgG1(E239D, M241L)
Reactivity	Human
Target	IGF-1R
Uniprot ID	P08069
Description	PE-conjugated Anti-IGF-1R(teprotumumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

