

PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 bM218 Clone ID

II 1884

CD218a; CDw218a; IL-1Rrp; IL1RRP; IL18RA Synonyme

Host Species Rabbit

PE-conjugated Anti-IL18RA antibody(DM218); Rabbit mAb Description

Delivery Under Development Uniprot ID IgG type Rabbit IgG Clonality Monoclonal Reactivity Human Flow Cyt Applications Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Liquid PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping Store at 2°C-8°C for 6 months

The protein encoded by this gene is a cytokine receptor that belongs to the interleukin 1 receptor family. This receptor specifically binds interleukin 18 (ILB); and is essential for ILB mediated singled transduction. This appear and ILB are reported to Induce the expression of this receptor in NR and TES. This gene along with four other members of the interleukin 1 receptor family; including ILBR; ILBR; ILBR; ILBR; ILBR (IL1R); and ILBR (T1:ST2); form a gene cluster on chromosome 2c, Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only Coniugate

PE-conjugated
All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

