

PRODUCT INFORMATION

Common Name	BI-655130, BI655130, BI 655130, spesolimab-sbzo
Synonyms	ILRL2; IL-36R; IL1RRP2; IL-36R
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid □ PBS with 0.05% Proclin300, 1% BSA
Host Species	Chimeric Humanized
IgG type	IgG1(K97R, L117A, L118A)
Reactivity	Human
Target	IL1RL2
Uniprot ID	Q9HB29
Description	PE-conjugated Anti-IL1RL2(spesolimab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

