

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 BMS-986015, IPH2102 |
| Conjugate | PE-conjugated |
| Synonyms | CD158b;CD158B1;NKAT-6;NKAT6;p58.2 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Homo sapiens |
| IgG type | IgG4 |
| Reactivity | Human |
| Target | KIR2DL2 |
| Uniprot ID | P43627 |
| Description | PE-conjugated Anti-KIR2DL2(lirilumab biosimilar) mAb |
| Delivery | Under Development |
| Yefei_Storage | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

DIMABIO CONFIDENTIAL

