

**PRODUCT INFORMATION**

<b>Clone ID</b>	5A5
<b>Target</b>	Mouse IgG
<b>Synonyms</b>	N/A
<b>Host Species</b>	Goat
<b>Description</b>	PE-conjugated Anti-Mouse IgG antibody(5A5), Goat mAb
<b>Delivery</b>	3-4 weeks
<b>Uniprot ID</b>	N/A
<b>IgG type</b>	Goat IgG
<b>Clonality</b>	Monoclonal
<b>Reactivity</b>	Mouse
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Purification</b>	Purified from cell culture supernatant by affinity chromatography
<b>Formulation &amp; Reconstitution</b>	Liquid□PBS with 0.05% Proclin300, 1% BSA
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	N/A
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated
<b>DIMA Disclaimer</b>	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

