

PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 bM88 Clone ID

NCR3;CD337;NKp30;1C7;LY117;MALS Synonyme

Host Species Rabbit

PE-conjugated Anti-NKp30 antibody(DM88); Rabbit mAb Description

Delivery Under Development Uniprot ID 014931 IgG type Rabbit IgG Clonality Monoclonal Reactivity Human Flow Cyt Applications Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Liquid PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping Store at 2°C-8°C for 6 months

The protein encoded by this gene is a natural cytotoxicity receptor (NCR) that may aid NK cells in the lysis of tumor cells. The encoded protein interacts with CD3-zeta (CD247); a T-cell receptor. A single nucleotide polymorphism in the 5' untranslated region of this gene has been associated with mild malaria suceptibility. Three transcript variants encoding different isoforms have been found for this gene.

Research use only Background

Conjugate PE-conjugated

.g of reverse All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

> Email: info@dimabio.com Website: www.dimabio.com



