

PRODUCT INFORMATION

| | |
|-----------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Common Name | Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 PF-04950615, RN-316, RN316 |
| Conjugate | PE-conjugated |
| Synonyms | FH3;FHCL3;HCHOLA3;LDLCQ1;NARC-1;NARC1;PC9 |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Humanized |
| IgG type | IgG2 |
| Reactivity | Human |
| Target | PCSK9 |
| Uniprot ID | Q8NBP7 |
| Description | PE-conjugated Anti-PCSK9(bococizumab biosimilar) mAb |
| Delivery | Under Development |
| Yefei_Storage | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

DIMABIO CONFIDENTIAL

