

PRODUCT INFORMATION

Common Name	IBI 306;IBI-306;IBI306
Synonyms	FH3;PC9;FHCL3;NARC1;LDLCQ1;NARC-1;HCHOLA3
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	N/A
IgG type	IgG2
Reactivity	Human
Target	PCSK9
Uniprot ID	Q8NBP7
Description	PE-conjugated Anti-PCSK9(tafolecimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

