

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\publicclass-woocommerce-print-products-public.php on line 2806 MEDI4736
Conjugate	PE-conjugated
Synonyms	PD-L1;CD274;B7-H1;PDCD1L1;PDCD1LG1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1µl/test
Formulation & Reconstitution	Liquid,PBS with,0.18%BSA,0.1%Proclin300
Host Species	Humanized
IgG type	IgG1(K97R, L117F, L118E, P214S)
Reactivity	Human
Target	PDL1
Uniprot ID	Q9NZQ7
Description	PE-conjugated Anti-PDL1 (durvalumab biosimilar) mAb
Delivery	In Stock
Yefei_Storage	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only




Figure 1. Flow cytometry analysis with 1 µl/test PE-conjugated Anti-PD-L1(durvalumab biosimilar) mAb on HEK293 cells transfected with Human PD-L1 protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

