

PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 bM87 Clone ID

PSCA;UNQ206;PRO232 Synonyme

Host Species Rabbit

PE-conjugated Anti-PSCA antibody(DM87), Rabbit mAb Description

Delivery Under Development 043653 Uniprot ID IgG type Rabbit IgG Clonality Monoclonal Reactivity Human Flow Cyt Applications Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Liquid PBS with 0.05% Proclin300, 1% BSA

Storage & Shipping Store at 2°C-8°C for 6 months

This gene encodes a glycosylphosphatidylinositol-anchored cell membrane plycoprotein. In addition to being highly expressed in the prostate it is also expressed in the bladder; placentar colory, idency, and stomach. This gene is up-regulated in a large proportion of prostate cancers and is also detected in cancers of the bladder and pancreas. This gene includes a polymorphism that results in an upstream start codon in some individuals; this polymorphism is thought to be associated with a risk for certain gastric and bladder cancers. Afternative splicing results in multiple transcript variants. Research use only

Email: info@dimabio.com Website: www.dimabio.com

Usage Coniugate

PE-conjugated
All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer



