

PRODUCT INFORMATION

Common Name	VX15/2503;VX 15;VX-15;VX 15;VX 15/2503;VX-15/2503
Synonyms	A8;GR3;BB18;CD100;COLL4;SEMAJ;coll-4;C9orf164;M-sema-G
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	N/A
IgG type	IgG4
Reactivity	Human
Target	SEMA4D
Uniprot ID	Q92854
Description	PE-conjugated Anti-SEMA4D(pepinemab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

