

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\publicclass-woocommerce-print-products-public.php on line 2806 XmAb-18087
Conjugate	PE-conjugated
Synonyms	SS-2-R;SS2-R;SS2R;SST2;SRIF-1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 2µl/test
Formulation & Reconstitution	Liquid,PBS with,0.18%BSA,0.1%Proclin300
Host Species	Chimeric Humanized
IgG type	IgG1
Reactivity	Human
Target	SSTR2
Uniprot ID	P30874
Description	PE-conjugated Anti-SSTR2(tidutamab biosimilar) mAb
Delivery	In Stock
Yefei_Storage	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only




Figure 1. Flow cytometry analysis with 0.25µl/test PE-conjugated Anti-SSTR2(tidutamab biosimilar) mAb on Expi293 cells transfected with Human SSTR2 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

