

PRODUCT INFORMATION

| | |
|---|--|
| Common Name | Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 D2E7 |
| Conjugate | PE-conjugated |
| Synonyms | DIF;TNF-alpha;TNFA;TNFSF2;cachexin;cachectin;TNFa |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Homo sapiens |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | TNFA |
| Uniprot ID | P01375 |
| Description | PE-conjugated Anti-TNFa(Adalimumab biosimilar) mAb |
| Delivery | Under Development |
| Yefei_Storage | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

DIMABIO CONFIDENTIAL

